Seminer

(Ücretsizdir)

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İ
Ç
İ
N
D
E
K
İ
L
E
R

Reprocessing of Medical Devices Cleaning / Disinfection / Sterilization			
Basics			
1. Legal aspects of sterilization and sterilization monitoring 2. Validation 3. Organization and documentation 4. EN ISO Standards and class 2 indicators according to EN ISO 11140-1 5. Scientific basics of process technology			
0.1-7	340		Dr. U. Kaiser 03/2011

Reprocessing of Medical Devices Cleaning / Disinfection / Sterilization			
Cleaning Processes Process Technology Routine Monitoring			
1. Critical variables influencing cleaning processes 2. Cleaning monitoring 3. gke Test method for cleaning indicators 4. Validation procedure of cleaning and disinfection processes according to EN ISO 15883 series 5. Standard EN ISO 15883 series for cleaning			
0.2-7	554		Dr. U. Kaiser 03/2011

Reprocessing of Medical Devices Cleaning / Disinfection / Sterilization			
Steam Sterilization Processes Process Technology Routine Monitoring			
1. Physical basics of the sterilization processes 2. Potential problems in steam sterilization process 3. Specific problems during sterilization of minimal invasive surgical (MIS-) instruments and tubes with small lumens 4. Routine monitoring 5. Use of medical device simulators (MDS) and batch monitoring systems (BMS)			
0.3-7	101		Dr. U. Kaiser 03/2011

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

MİC* uçları

İkiz uçlar

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- » K- Tel ve dispenseri
 - Tek-Yuvalı
 - » K-Teli depola, dağıt ve and sterilize et
 - » Tek seferde dağıt
 - » Perforated for proper sterilization
 - » Clearly marked with diameter for quick identification
 - » Two lengths: 6.25” and 12.25”
 - Çoklu-Yuvalı
 - » Store, sterilize and dispense four different sizes
 - » Dial allows for easy dispensing one at a time
 - » Chambers are clearly marked for easy identification

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Legal aspects for the implementation of the European Medical Device Directive (MDD) in connection with validation, monitoring and documentation of sterilization processes

In all countries of the European Community and in the associated member countries the Medical Device Directive (MDD) is finally implemented since 1998. All sterile products being used at a human being are defined as medical devices according to the MDD. To produce sterile products the validation of a sterilization process is required. It is not sufficient to commission the sterilizer - often misinterpreted as "validation". A sterilization process consists of the following components:

- (1) sterilizer with all accessories and supply ingredients
- (2) sterilization program selected
- (3) type of goods sterilized
- (4) packaging material
- (5) load configuration

If one of those components is modified (especially other types of goods) a new validation of the process is required.

The validation of a sterilization process is a procedure (see point 1 below) generating a test report which guarantees that all goods being sterilized have achieved a sterility assurance level (SAL) of 10^{-6} CFU per part.

If those goods are sold to a third party (brought to the market) according to the MDD the validation must be checked by a Notified Body. This process is called "conformity process". According to the second modification of the MDD health care facilities that re-process their sterile goods do not bring products to the markets. The re-processing is classified as a "maintenance process". Therefore an individual conformity process by a notified body is not required in this case.

However, for the application of those goods used at the patient, health care facilities require the same validation, monitoring and documentation procedures which have to be carried out under the responsibility of the health care facility:

(1) validation

The validation has to be carried out as follows with the steps below:

1. Installation and operation qualification (IQ and OQ) called commissioning of the sterilizer with a type test according to the specification of the manufacturer
2. Performance qualification (PQ) checking the sterilization process including all goods, packaging etc. by parameter testing, the use of process challenge devices (PCDs) with chemical and/ or biological indicators to proof the sterilization in all worst case areas of the load

Under certain sterilization conditions for complex goods, parametric release may not be sufficient and biological indicator suspensions have to be used at worst case places for monitoring and validation.
3. Validation report of the commissioning and PQ according to the technical state of the art.

Sterilization as part of an industrial production process, i.e. production of sterile disposables, can be easily validated because the sterile goods are

and therefore the whole sterilization process is always identical.

Sterilization processes in health care, i.e. in hospitals, are more complex because the goods to be sterilized and the load configuration always differ. The validation of this process therefore is much more difficult. In order to prove the success of the sterilization process the validation has to refer to the so-called worst-case load configuration. Concerning steam sterilization the validation standard EN ISO 17665-1 contains the requirement to define the worst-case load (that means the configuration that is most difficult to sterilize) before any measurement starts. Under the condition all other load configurations are easier to sterilize, the worst-case load can be used for validation in health care units.

(2) Routine monitoring

Routine monitoring has to be carried out during every cycle to proof that each cycle is running according to the validated process. Typically temperature and pressure-time-diagrams and also the penetration of the sterilant have to be monitored with a representative process challenge device (PCD).

(3) Documentation

Most parameters of a process can be monitored by measuring the physical data. Any manual part of a process has to be standardized in order to ensure that it always is identical and therefore comparable.

A wide-spread means of standardization is provided by the European standard EN ISO 9001. Installing a quality management system can be used to organize the standardized processes, i.e. in so-called SOPs (Standard Operation Procedures).

All relevant data collected during the monitoring process have to be documented at a documentation chart i.e. the batch number for each cycle, the production date, the responsible person, and – if applicable – the expiry date. The same data but at least the batch number has to be documented on each pack.

In case a patient gets a secondary infection, legal proceedings may be raised and the health care facility has to provide data to proof that the re-processing of sterile goods has been carried out according to the latest technology standard.

This standard is reflected in local laws and in European and international standards like:

EN-ISO 14937

Validation and routine monitoring for sterilization processes

EN ISO 11137-1

Validation of radiation sterilization processes

EN 15424

Validation of Low-Temperature Steam Formaldehyde (LTSF) sterilization processes

EN ISO 11135-1

Validation of ETO sterilization processes

EN ISO 17665-1


Validation of steam sterilization processes

EN 556-1 + -2

Definition of sterility assurance level

The *gke* application laboratory may provide more detailed information about validation, monitoring and documentation of your individual sterilization process.

Processing Cycle of reusable sterile Goods

9.1-4

259

U. Kaiser

07/2005

Requirements for a Quality Management System in the CSSD-Department (Extract of ISO 9000:2000 and ISO 13485)

- Written standard operation procedure for all relevant jobs (SOPs)
- Assignment of all SOPs to the responsible person in charge (job description)
- Substitution plan
(Also the substitute person requires a job description.)
- Training of all responsible people in charge (also the substitute persons)
- Monitoring of all relevant process parameters
- Batch related documentation of relevant process
- Patient related documentation based on the batch produced

9.2-4

260

U. Kaiser

07/2005

General information about National, European and International standards

National Standards

DIN - Deutsche Industrie Norm
BS - British Standard etc.

European Standards

are mandatory for all EU member countries including the associative countries Norway, Iceland and Switzerland. EU member countries are: Austria, Belgium, Bulgaria, Czech Republic, Cyprus, Denmark, Estonia, Finland, France, Germany, Greece, Hungary, Ireland, Italy, Latvia, Lithuania, Luxembourg, Malta, Netherlands, Poland, Portugal, Romania, Slovakia, Slovenia, Spain, Sweden, and the United Kingdom.

CEN - Comité Européen de Normalisation
EN XXX - European Norm
prEN XXX - Preliminary European Norm (Draft, no final Standard)

In case a European standard is published, national standards with the same content automatically become obsolete within the EC member countries.

World Wide Standards

ISO - International Standard Organization
ISO XXX - International Standard
ISO/DIS XXX - Draft International Standard (no final version)

ISO members are the national standard organizations of all European countries, USA, Japan and many other countries.

In 2005 many ISO- and EN standards merged to one identical EN-ISO standard.

Standards for biological and chemical indicators as well as standards for validation are currently merging.

EN-ISO XXX - International Standard, taken over as European Norm as well

10.1-29

263

U. Kaiser

03/2008

Standards for Reprocessing Processes (1)

Validation	Sterilizers	Washing Disinfection	Chemical Indicators	Biological Indicators	Packaging
EN ISO 14937 Requirements for development, validation and routine monitoring of all sterilization processes	EN 285 Requirements for large sterilizers (over 54 l)	EN ISO 15883-1 General Requirements for washer/disinfectors	EN 867-5 Chemical indicator systems for steam sterilizers (Test standard for hollow load test)	EN ISO 11138-1 General requirements and classifications on Biological Indicators (BI)	EN ISO 11607-1 Packaging of medical devices
EN ISO 11135-1 -2 EO processes	EN 13060 Requirements for small sterilizers (below 54 l)	EN ISO 15883-2 W/D requirements for surgical instruments	EN ISO 11140-1 General requirements, classifications and test procedure for chemical indicators (CI)	EN ISO 11138-2 BI for EO sterilization	EN ISO 11607-2 Validation requirements for forming processes
EN ISO 11137-1 -3 Radiation processes	EN 14180 Requirements for LTSF sterilizers	EN ISO 15883-3 W/D requirements for containers for human waste	EN ISO 11140-3 Requirements for the original BD-test page	EN ISO 11138-3 BI for steam sterilization	prEN ISO/PTS 16775 Guidance for the application of EN ISO 11607-1+2
ISO/PDTS 13004 Radiation processes	EN 1422 Requirements for EO sterilizers	EN ISO 15883-4 W/D requirements for thermolabile endoscopes	EN ISO 11140-4 Requirements for BD-Simulation tests	EN ISO 11138-4 BI for dry heat sterilization	EN 868 Series 2-10 Packaging of sterile goods
EN ISO 17665-1 -3 Steam processes	EN ISO 18472 Requirements for test sterilizers (resistometers)	EN ISO 15883-5 W/Ds - test soils and methods	ISO 11140-5 Requirements for the US BD-test	EN ISO 11138-5 BI for LTSF sterilization	
EN 15424 = EN ISO 25424 LTSF processes	DIN 58951 Requirements for steam sterilizers in laboratories	EN ISO 15883-6 W/Ds - Requirements and tests for general purpose W/Ds with thermal disinfection	prEN ISO 11140-6 will substitute EN 867-5 Requirements for Hollow load test	ISO 11138-6 BI for H ₂ O ₂ sterilization (under development, no draft available)	
EN ISO 14937 also for H ₂ O ₂ / Plasma processes, since no special standard available	DIN 58946-7 Requirements for installation and operation of steam sterilizers	prEN ISO 15883-7 W/Ds - Requirements and tests for general purpose W/Ds with chemical disinfection for bedframes, containers, etc.	EN ISO 15882 Guidance for the selection, use and interpretation of the results for chemical indicators	EN ISO 14161 Guidance for the selection, use and interpretation of the results for biological indicators	
prEN ISO 20857 Dry heat processes	DIN 58948-7 Requirements for installation and operation of low temperature sterilizers				
EN ISO 17664 Information about reprocessing of re-usable medical devices					
DIN 58921 Validation of medical device simulators (MDS) (English version available)					
EN 556-1 Definition: Sterility Assurance Level					

European Medical Device Directive (MDD) 93/42

10.2-29

264-1

U. Kaiser

02/2013

Standards for Reprocessing Processes (2)

Pharmaceutical Procedures	Sterilizing agents	Disinfectants and disinfectors	Aseptical Production	Additional standards	
DIN 58950-1 Definitions	EN ISO 14160 Liquid chemical sterilizing agents for medical devices	EN 1499 Hygienic cleaning of hands	EN ISO 13408-1 General Requirements	EN 980 Symbols for labeling of medical devices	EN ISO 11139 Terms and definitions in sterilization standards
DIN 58950-2 Technical requirements		EN 1500 Hygienic hand disinfection	EN ISO 13408-2 Filtration	EN 1041 Information supplied by the manufacturer of medical device	EN ISO 11737-1 - 2 Microbiological methods
DIN 58950-3 Tests		DIN 12353 Preservation of test organisms	EN ISO 13408-3 Lyophilization	EN 15224 Healthcare services	EN ISO 14971 Risk management of medical devices
DIN 58950-6 Operation		EN 14476 Virucidal quantitative suspension test for chemical disinfectants and antiseptics	EN ISO 13408-4 Clean-in-place technologies	EN ISO 13485 Medical device quality management system	EN 15986 Symbols to mark medical devices
DIN 58950-7 Requirements on services and local environment		DIN 59949 Steam disinfection apparatus	EN ISO 13408-5 Sterilization in place	EN 15223-1 Symbols for labeling of medical devices	DIN 58953 Sterile supply – Terminology, Logistics
		RKI ¹ list of tested disinfectants and disinfection processes	EN ISO 13408-6 Isolator systems	EN ISO 10993-1 -17 Classification of medical devices	
	Surgical drapes	VAH ² list of disinfectants			
	EN 13795 Surgical drapes, gowns and clean air suits General requirements				

¹ RKI = Robert Koch Institute, Germany

² VAH = Association for applied hygiene, Germany

10.3-29

264-2

U. Kaiser

02/2013

Merged European CEN- and international ISO-standards for non-biological (chemical) indicators

European Standards EN 867, parts 1 - 5 Non-biological systems (chemical indicators) for testing sterilization processes. (will be withdrawn)		ISO-Standards 11140:1994, parts 1 - 2: Non-biological systems (chemical indicators), (old version is withdrawn)		(EN) ISO 11140:2005, parts 1 - 6 New standards combining old EN and ISO standards	
Part 1	General requirements	Part 1 withdrawn	Classification and general requirements class 1-6	Part 1	Classification and general requirements class 1-6
Part 2 withdrawn	Process indicators (class A) for use in steam, EO, FO and dry heat	Part 2 withdrawn	Requirement on test equipment for steam, EO, FO and dry heat sterilization processes	Part 2 does not exist	EN-ISO 18472 test equipment for biological and chemical indicators (separate Standard)
Part 3 withdrawn	Class B indicators-sheets for use in the Bowie-Dick-Test	This Euro norm is transferred to EN-ISO Standard and is withdrawn as EN Standard		Part 3	Indicators of class 2 for test sheets for steam penetration (= EN 867-3)
Part 4 withdrawn	Class B indicators for use as alternative to the BD-Test	This Euro norm is transferred to EN-ISO Standard and is withdrawn as EN Standard		Part 4	Indicators of class 2 for test packs for steam penetration; reference BD-test pack according to the European standard EN 285 (= EN 867-4)
New, no EN- or ISO Standard existing until now. The new standard is no EN but only an ISO standard				Part 5	Indicators of class 2 for test sheets and test packs for the air removal; test reference BD-test-pack according to AAMI – USA
Part 5 will be withdrawn	Class B indicator system for use in large sterilizers according to EN 285 and small sterilizers according to EN 13060	This Euro norm is currently modified into an EN-ISO Standard, after that EN 867-1 and -5 are withdrawn completely		Part 6 (Draft)	Class 2 indicator for use in large sterilizers acc. to EN 285 and table- top sterilizers acc. to EN 13060 (= EN 867-5)
Will be revised and adopted to the new EN-ISO 11140 series				<u>EN ISO 15882 Chemical indicators</u> Guidance for selection, use and interpretation of results	

10.5-29

266

U. Kaiser

07/2009

European Norms (EN) and/or International Standards (ISO) for sterilizers, packaging material, sterilization monitoring and validation

7. Technical requirements for sterilizers

- EN 285:2008 Healthcare steam sterilizers > 54 l
To the BD type test the hollow load type test according EN 867-5 had been added to assure that hollow devices and MIS instruments can be sterilized
- EN 13060 Small steam sterilizers < 54 l

EN 13060 divides small steam sterilizers into the following types:

Steam sterilization process	Application	Type
Fractionated vacuum	all types of packaged, hollow, solid and porous goods have to be tested with the hollow load PCD	B
Gravity displacement or single vacuum	non-packaged, solid goods to be used immediately after sterilization or for goods to be sterilized to prevent cross-infections (no hollow devices)	N
Depending on the manufacturer's specification	goods, specified by the manufacturer	S

Type B sterilizers require fractionated vacuum and must pass the hollow load test (EN 867-5).

- EN 14180 Low temperature steam-formaldehyde sterilizers
Requirements and testing methods
- EN 1422 Ethylene oxide sterilizers – Requirements and testing methods
- EN ISO 18472 Resistometer (Test sterilizer) for biological and chemical indicators
- DIN 58951 Steam sterilizers for laboratories - Part II: Requirements
- All sterilizers (re-)sterilizing medical products are medical devices (MDs) and classified as 2b in the Medical-Device-Directive (MDD).

10.12-29

271

U. Kaiser

07/2009

Differences between small and big steam sterilizers for the health care market

Usable volume 30x30x60 cm = 1 STU	Sterilizer standard	Requirement of the standard / type test
< 54 l or > 54 l 1 STU does not fit into the chamber	small sterilizers EN 13060 type - B - N - S	Successful air-removal hollow load helix test according to EN 867-5 According to the specifications of the manufacturer
≥ 54 l 1 STU or more fits into the chamber	big sterilizers EN 285	<u>Edition 1994:</u> BD-Test 7 kg (European) <u>Since 01/2008:</u> BD-Test 7 kg, <u>and</u> hollow load test according to EN 867-5

10.13-29

272

U. Kaiser

07/2009

Classification of chemical indicators for different sterilization processes according to the International standard EN-ISO 11140-1

Class	Description
I	Process indicators
II	Indicators for specific use (Details are found in other Standards)
III	Single parameter indicators
IV	Multi-parameter indicators
V	Integrating indicators (emulating biological indicators)
VI	Emulating indicators

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Quote in attachment of annex 4
„Special information for the use of chemical indicators“
(page 1273-1274)

Page 1 of EN ISO 11140 differentiates several classes of chemical indicators with different challenges. However this classification is no hierarchical graduation („better“ or „worse“).

[...]

Indicator systems of class 2 are required for special tests, which are required in the standards about sterilizers, they consist of a Process Challenge Device (PCD) and an indicator inside.

10.16-29

369

J. Metzinger

02/2013

Definition and construction of Indicators according EN ISO 11140-1

Indicators of the classes 1, 3, 4, 5 and 6:

- Consist of indicator substance on a carrier material (eventually self-adhesive indicator strip)
- All Specifications are determined in EN ISO 11140-1 for the different indicator classes and sterilization processes

10.17-29

275

J. Metzinger

04/2008

Definition and construction of Indicators according EN ISO 11140-1

Indicators of class 2 according EN ISO 11140-1:

- Consist of a specific test load (PCD = Process Challenge Device) and indicator system (indicator strip) or another detector
- The specifications of class 2 Indicators are not described in EN ISO 11140-1, but in other standards, i.e.:

EN 285 (Bowie-Dick-Test),
EN ISO 11140-3 (Indicator system / test sheet for Bowie-Dick-Test in EN 285,
EN ISO 11140-4 (test method to validate BD-Simulation Tests),
ISO 11140-5 (USA BD-Simulation Test)
EN 867-5 (hollow load helix test) - standard currently modified to become
Draft EN ISO 11140-6
DIN 58921 (Medical Device Simulator)
etc.
- A class 2 indicator is a combination of test device (PCD) and indicator system. The single components cannot be used separately without revalidation and are not defined in a class (see picture).
The question: „Which class of indicator is inside the PCD?“
is an invalid question.

10.18-29


276

J. Metzinger

11/2010

Definition and construction of the indicator according EN-ISO 11140-1

10.19-29

277

gke

J. Metzger

03/2008

Difference between habitual language use and definitions according EN-ISO 11140-1 class 2 indicators

Habitual language use	Picture	Name according EN-ISO 11140-1
Indicator strip to be used in PCD		Indicator system
Process Challenge Device (PCD) (without indicator strip)		Specific test load
Indicator system (PCD with indicator strip inside)		Indicator (class 2 Indicator) in the form it is used

10.20-29


278

gke

J. Metzger

07/2007

**PTFE-Helix-PCD design according
EN 867-5 and EN ISO/DIS 11140-6**

10.23-29

333

gke

U. Kaiser


12/2010

**Specification and handling of
class 2 indicators according EN ISO 11140-1**

c) Requirements on the leak tightness

Seal test according EN 867-5, paragraph 5.1.6.1
Seal of demountable cap on process challenge device

The seal test describes as test setup a PCD which builds up pressure inside using a syringe and an adapter. The closed PCD is put into a water bath (according 5.1.6.1) and into an oil bath (according 5.1.6.2) with maximum sterilization temperature, i.e. 140°C. When pressurized with a syringe, no bubbles must be visible.

10.26-29

283

gke

J. Metzger

04/2008

Pathogenic micro-organisms causing disease

1. Micro-organisms (with DNA) and Exotoxins

Bacteria	(micro-organisms)
Protozoa	(animals)
Fungi	(plants)
Viruses	(obligate cell parasites of bacteria, animals, plants)
Exotoxins	toxic substances that are released from bacteria (e.g. enterotoxins, diphtheria toxins, botulinus toxins)

These substances can be inactivated in standard sterilization processes.

2. Chemical substances (without DNA)

Prions	Proteins responsible for BSE, CJD (Creutzfeld Jacob disease)
Pyrogens	Endotoxins toxic substances inside of bacteria which are released during decomposition of mainly Gram-negative bacteria (e.g. Lipopolysaccharides LPS) Chemical substances, that cause an immune defense reaction

Those substances cannot be securely inactivated in standard sterilization processes.

Resistance of important micro-organism-groups depending on temperature and time

	Temperature of water	80°C	100°C	121°C	134°C
Ia	Plasmodia Flagellats Viruses Bacteria without spores Yeasts Moulds	1 – 5 min	--	--	--
Ib	Spores of yeasts and moulds	5 – 10 min	1 min	--	--
II	Spores of germs with low resistance	--	1 – 60 min	1 min	--
III Reference indicator	Spores of germs with higher resistance (<i>G. stearothermophilus</i>)	--	60 min – 60 hrs.	8 min	1 min
IV	Thermo-resistant spores	--	--	--	up to 6 hrs.
V	Priones (Kreutzfeld-Jacob-disease)				30 – 60 min

Standardized Biological Organisms for Sterilization monitoring

Bacteria type	Sterilization process to be checked	ATCC-No.: *	Incubation temperature [°C]
1. <i>B. atrophaeus</i> (formerly : <i>B. subtilis</i>)	Dry heat, Ethylene oxide, Formaldehyde Steam } for disinfection only	9372	35°C
2. <i>G. stearothermophilus</i>	Steam, Formaldehyde Hydrogen peroxide	7953	55°C
3. <i>B. pumilus</i>	γ - and cobalt radiation	27142	35°C

* American-Type-Culture-Collection

2.3-17


110

U. Kaiser

07/2005

The survival curve of microorganisms per part for the steam sterilization process at different temperatures

2.4-17


111

U. Kaiser

07/2005

G. Stearotherophilus Temperature - Time - Dependent - Kill - Kinetics

2.5-17

112

gke

U. Kaiser

09/2007

Sterility Assurance Level (SAL)

$$\frac{1}{10} = \frac{\text{Number of contaminated packs}}{\text{Number of total packs}} = 0,1 = 10^{-1} [\text{CFU/part}]$$

$$\frac{10}{100} = \frac{\text{Number of contaminated packs}}{\text{Number of total packs}} = 0,1 = 10^{-1} [\text{CFU/part}]$$

$$\frac{1}{100} = \frac{\text{Number of contaminated packs}}{\text{Number of total packs}} = 0,01 = 10^{-2} [\text{CFU/part}]$$

•
•
•

$$\frac{1}{1.000.000} = \frac{\text{Number of contaminated packs}}{\text{Number of total packs}} = 0,000001 = 10^{-6} [\text{CFU/part}]$$

S
A
L

}

terility
assurance
level

EN 556-1: products sterilized in packs: SAL=10⁻⁶ [CFU/part]

EN 556-2: sterile filling of liquids: SAL ≥ 10⁻³ [CFU/part]

2.6-17

113

gke

U. Kaiser

06/2008

Total resistance of biological indicators

- is depending on: **(1)** logarithm of the population (lgPop)
 (2) resistance of the germs (D-value)

Total resistance of a BI = F_{BIO} -value [min] = lgPop x D-value [min]

The F_{BIO} -value is the time to reduce the population
to the mean value of 1 CFU (66% growth)

% Growth	Sterilization time	Population [CFU]	D-value [min]	F_{BIO} [min]
100	(lgPop-2) x D-value = survival	10 ⁶	1.5	9
90	(lgPop-1) x D-value	10 ⁵	2.0	10
66	lgPop x D-value			
10	(lgPop+1) x D-value			
1	(lgPop+2) x D-value			
0.1	(lgPop+3) x D-value			
0.01	(lgPop+4) x D-value = kill			

The indicator with the lower population has a
higher total resistance but is less expensive.

2.7-17


114

U. Kaiser

06/2011

Determination of the z-value

The z-value (expressed in °C) is the temperature difference which is necessary to change the D-value by factor 10.

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
05/2012

Colour Reference Chart

for all **gke Clean-Record®** Cleaning indicators

Actions for troubleshooting

Worse test results, compared to earlier batches, can occur in different cases:

Possible causes	Action
Wrong program	Check documentation! Check if the correct program for the load has been used?
Modified load configuration	Check, if the load configuration complies with the specifications documented in the quality management system and the validation report.
Place of cleaning indicator	Check, if the cleaning indicator has been placed at the right location.
Other cleaning detergent	Check cleaning detergent can. Did you use the right detergent and is it expired?
Wrong dosage	Check dosage mark. Mark liquid level in cleaning detergent can, run the program and check the liquid level or weight difference!
Incorrect temperature curve	Compare temperature curve from current batch with previous ones.
Spray arm not moving	Check, if the spray arm can be turned without resistance. Machine with glass door: Check spray arm during cleaning process. Other machines: Stop program during process and check if spray arm is able to move.
Dirt strainer is clogged	Check dirt strainer, clean if required.
Modified water quality	If tap water is used: Check hardness and salt content. Check if values have changed. If softened or demineralized water is used: Check water softening or demineralization system, pH-value and conductivity.

After checking the above mentioned actions, repeat cleaning process.
Modification of spray and flush dynamics, e.g. clogged nozzles, leakages etc. may not be detected by yourself. To detect these possible causes, call technical service.

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Cleaning Agent	Alkaline detergent, 0,5 %, pH = 7,7, with enzymes					
Water quality	Demineralized water					
Flow rate	1,0 l/min			3,0 l/min		
Split width HF-PCD	2 mm	3 mm	4 mm	2 mm	3 mm	4 mm
Flow time						
10 sec						
30 sec						
60 sec						
3 min						
5 min						
10 min						

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02/2013

gke
Cleaning and Sterilization Monitoring

Hemp/Clasp (<i>Showering Hospital</i>)		Department: <i>6250</i>	City: <i>Excellist II</i>
Manufacturer: <i>Paul Chen Ltd</i>	Unit No: <i>1276</i>		
Program No: <i>007</i>	Batch No: <i>1</i>	Date: <i>2012-04-26</i>	
Fixation of Indicator: <i>gkx</i> ☒ in the below flow PCD ☒ 2 mm Split width (green) ☒ 3 mm Split width (blue) ☒ 4 mm Split width (red) ☐ other location:	Cleaning description: Manufacturer: <i>Egypt Chemical</i> Product Name: Best before: <i>12/2013</i>	Adhesive cleaning indicator here:  Batch No. of Indicator: <i>1206 1144</i>	
Approval: ☒ yes ☐ no	Checked by: <i>Paul Green</i>		
Program No: <i>104</i>	Batch No: <i>2</i>	Date: <i>2012-04-26</i>	
Fixation of Indicator: <i>gkx</i> ☒ in the below flow PCD ☒ 2 mm Split width (green) ☒ 3 mm Split width (blue) ☒ 4 mm Split width (red) ☐ other location:	Cleaning description: ☒ as above Manufacturer: Product Name: Best before:	Adhesive cleaning indicator here:  Batch No. of Indicator: <i>1270 1140</i>	
Approval: ☒ yes ☐ no	Checked by: <i>Paul Green</i>		
Program No: <i>052</i>	Batch No: <i>3</i>	Date: <i>2012-04-20</i>	
Fixation of Indicator: <i>gkx</i> ☒ in the below flow PCD ☒ 2 mm Split width (green) ☒ 3 mm Split width (blue) ☒ 4 mm Split width (red) ☐ other location:	Cleaning description: ☒ as above Manufacturer: Product Name: Best before:	Adhesive cleaning indicator here:  Batch No. of Indicator: <i>1260 1200</i>	
Approval: ☒ yes ☐ no	Checked by: <i>Paul Green</i>		
Program No: <i>003</i>	Batch No: <i>4</i>	Date: <i>2012-04-26</i>	
Fixation of Indicator: <i>gkx</i> ☒ in the below flow PCD ☒ 2 mm Split width (green) ☒ 3 mm Split width (blue) ☒ 4 mm Split width (red) ☐ other location:	Cleaning description: ☒ as above Manufacturer: Product Name: Best before:	Adhesive cleaning indicator here:  Batch No. of Indicator: <i>1260 1200</i>	
Approval: ☒ yes ☐ no	Checked by: <i>Paul Green</i>		

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**Reprocessing of Medical Devices
Cleaning / Disinfection / Sterilization**

**Steam Sterilization Processes
Process Technology
Routine Monitoring**

1. Physical basics of the sterilization processes
2. Potential problems in steam sterilization process
3. Specific problems during sterilization of minimal invasive surgical (MIS-) instruments and tubes with small lumens
4. Routine monitoring
5. Use of medical device simulators (MDS) and batch monitoring systems (BMS)

Overview of common used sterilization processes

Thermal		Chemical		Physical	
Process	Temperature [°C]	Process	Temperature [°C]	Process	Temperature [°C]
Dry heat	160 – 180	Ethylene oxide	30 – 70	• γ radiation sterilization	20
Steam with air removal: • gravity displacement • single vacuum • vacuum-steam-injection • high-vac • fractionated vacuum	110 – 135	Formaldehyde	40 - 70	• β radiation sterilization	20
		Hydrogen peroxide (Plasma)			30 – 50
Microwave	110 - 135				

Definition "Sterile":

Complete inactivation of viable germs on a product


**The number of alive germs can be reduced continuously with increasing sterilization time but SAL cannot reach zero.
(Theory of kill kinetics)**

According EN 556-1 in the EU a product can be labeled sterile, if the sterility assurance level (SAL) is $\leq 10^{-6}$.


EN 556-2 accepts $SAL \leq 10^{-3}$ when sterile liquid filling is done.

In the USA the required SAL is depending on the application risk ($SAL \leq 10^{-3} - 10^{-9}$).

Thermal Sterilization – Dry Heat Sterilizer

Principles of the gravity displacement sterilization cycle (1)

3.3-55


125

gke

U. Kaiser

07/2005

Principles of the gravity displacement sterilization cycle (2)



3.4-55

126


gke

U. Kaiser

07/2010

Trans-atmospheric Fractionated Vacuum Sterilization Process

3.15-55

133

gke

U. Kaiser

11/2005

									No.	Goods to be sterilized	Packing	Difficulty of Air-Removal
+ successful sterilization process									1	Instruments without holes	without	easy
- unsuccessful sterilization process									2	Instruments without holes	with	
									3	Porous goods (Cotton etc.)	with or without	
									4	Hollow instruments like MIS*-instruments, endoscopes, catheters etc.	with or without	difficult
No.	Sterilization Process	Type of small sterilizers	Process-Characteristic	Air-removal					Sterilization-Process-Quality-Dependence from the Steam-Sterilization-Process and the Goods sterilized *MIS = Minimal-Invasive-Surgical			
1	Gravity-Cycle (Gravity-Displacement)	N	Displacement	poor	+	-	-	-				
2	Single Vacuum Cycle 100 mbar	S			+	(+)	-	-				
3	Vacuum-Injection-Cycle	S	Steam		+	+	+	-				
4	Super-Atmospheric Pressure Cycle (Over Pressure Cycle)	S			+	+	+	-				
5	Vacuum-Over-Pressure Cycle	S			+	+	+	-				
6	Vacuum-Over-Pressure Cycle (stairs)	S			+	+	+	-				
7	Sub- or trans-atmospheric Pressure Cycle (Fractionated Vacuum)	B		good	+	+	+	+				

3.19-55

137

gke

U. Kaiser


01/2006

Which are the necessary sterilization parameters in steam sterilization processes? (1)

- 1 - Steam
- 2 - Condensation of Steam to Water
- 3 - Pressure
- 4 - Temperature-Time-Window

Dry Cellulose Packs (Cotton or Paper) create superheated non-condensing steam by hydroscopic condensation

Biological indicators are not inactivated under superheated steam conditions.

**Which are the necessary sterilization parameters
in steam sterilization processes? (2)**

1 - Steam?

Steam is water in the gaseous phase. Germs and biological indicators are not killed at the same speed in steam superheated steam using the same temperature-time-window (i.e. 121°C - 15 min; 134°C - 3 min) in comparison to the speed in steam sterilization processes.

The kill velocity in non-condensed steam is similar to dry heat or in anhydrous liquids and requires similar temperature-time-windows

(i.e. 160°C - 2 h; 180°C - 0.5 h).


- NO -

**Steam Sterilization of infusion solution condensing steam at
the outside glass walls**

(no condensation inside the bottle)

2 - Condensation of Steam to Water?

Alternatively microwaves, infrared radiation, heating plates can be used to heat up the water

**Which are the necessary sterilization parameters
in steam sterilization processes? (3)**

2 - Condensation of Steam to Water?

Condensing steam creates an excellent heat transfer to all areas where steam can condense to water. However, liquid sterilization in water demonstrates that the kill velocity of germs is identical in water without any condensation under the same temperature-time-window like in steam sterilization processes.

During the plateau period in steam sterilization processes there is no condensation, however the same kill velocity occurs like in water.

- NO -

3.24-55

142

U. Kaiser

01/2011

**Which are the necessary sterilization parameters
in steam sterilization processes? (4)**

3 - Pressure?

Pressure in the scale region of 1 – 10 bar or 1.000 kPa does not influence the kill velocity of germs in any sterilization process.

Pressure in steam sterilization processes is required to achieve the necessary temperature. The boiling point of water depends on the pressure (see table below).

Pressure [bar]	Temperature [°C]
1	100
2	121
3	134

- NO -

3.25-55

143

U. Kaiser

01/2011

**Which are the necessary sterilization parameters
in steam sterilization processes? (5)**

4 - Temperature-Time-Window?

If the temperature-time-window used in steam sterilization processes (121°C - 15 min; 134°C - 3 min) is also used in dry heat or anhydrous liquids, the kill velocity of germs or biological indicators is much lower than in steam sterilization processes under the same temperature-time-window.

- NO -

3.26-55

144

U. Kaiser

01/2011

**Temperature and time parameters to achieve overkill for the
steam sterilization process according to EN ISO 17665-1**

4 - Temperature-Time-Window?

Temperature [°C]	Sterilization time [min]*	Equilibration time [min]**	F _{0 121°C} [min]	Remarks
121	15	<0,5	15	Those conditions are only valid in presence of water (wet surfaces).
134	3	<0,5	ca. 60	

**Dry heat and superheated steam sterilization process
(incl. non-polar solvents and oils)**

Temperature [°C]	Sterilization time [min]*	Equilibration time [min]	Remarks
160	120	10 – 50	The temperature-come-up-time changes depending on the heat capacity of the goods and the insulation of the pack.
180	30	10 - 30	

* Sterilization time after reaching the temperature of the goods on the surface and in hollows

** If complete air removal is achieved

3.27-55

145

U. Kaiser

05/2011

Which are the necessary sterilization parameters in steam sterilization processes? (6)

In liquid sterilization with pure water the kill velocity of germs or biological indicators is identical like in steam sterilization processes using the same temperature-time-window. The resistance testing of biological indicators for steam sterilization processes is carried out in water to prevent super heat conditions.

The sterilization agent in steam sterilization processes is:

- **Water**, having a defined temperature-time-window applied to the surface of the instruments (e.g. 121°C – 15 min; 134°C – 3 min = F₀-value)
- All surfaces which should become sterile, have to be wet

This information is written already since 10 years in the indicator standard EN ISO 11140-1 in 5.2:

For the different sterilization processes, the following variables are defined as being critical:

Steam	Time, temperature and water (as delivered by saturated steam)
Dry Heat	Time and temperature
EtO	Time, temperature, relative humidity and ethylene oxide (EO) concentration
Steam-FO	Time, temperature, water (as delivered by saturated steam) and formaldehyde concentration
Vaporized H ₂ O ₂	Time, temperature, hydrogen peroxide concentration, and, if applicable, plasma

3.28-55

146

U. Kaiser

01/2011

Reasons for superheated steam

Superheated steam does not contain any condensate in the gas phase and is present in a steam sterilizer under the following conditions:

1. Pressure reduction in the steam pipe without sufficient cooling before entering the sterilizer
2. Sterilizer walls heated up to a higher temperature than the saturated steam temperature inside the sterilizer
3. Hygroscopic condensation inside non-conditioned cellulose materials (like cotton balls, packaging material, cotton and linen)

Steam in the gasphase does not sterilize at a temperature up to 134°C, only water can sterilize under those temperature and time conditions.

3.29-55

359

U. Kaiser

09/2012

Steam consumption and contraction during condensation

about 350 – 400 l/10 kg load

Steam

1000 ml steam → approx. 1 ml condensate (water)

Significant volume contraction of steam due to condensation

- The heat transfer occurs directly at the condensation location on the instruments
- During condensation the steam volume diminishes (converted to water), consequently additional steam is forced into the generated vacuum

3.30-55


147

gke

U. Kaiser

01/2011

Volume contraction of water during Condensation Non-condensable gases (NCG) do not reduce their volume

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
U. Kaiser

01/2011

Temperature and time parameters to achieve overkill for the steam sterilization process according to EN ISO 17665-1

Temperature [°C]	Sterilization time [min]*	Equilibration time [min]**	F ₀ 121°C [min]	Remarks
121	15	<0,5	15	Those conditions are only valid in presence of water (wet surfaces).
134	3	<0,5	ca. 60	


Dry heat and superheated steam sterilization process (incl. non-polar solvents and oils)

Temperature [°C]	Sterilization time [min]*	Equilibration time [min]	Remarks
160	120	10 – 50	The temperature-come-up-time changes depending on the heat capacity of the goods and the insulation of the pack.
180	30	10 - 30	

* Sterilization time after reaching the temperature of the goods on the surface and in hollows

** If complete air removal is achieved

Comparison of separations of non-condensable gases (NCG) in porous loads and hollow instruments

3.36-55


153

gke

U. Kaiser

07/2005

Internal volume in tubes and minimal invasive surgical (MIS-) instruments Thermo-sensor testing inside of metal tubes

- (1) Non-condensable gases (NCG) < 1 ml become already critical in minimal invasive surgical (MIS-) instruments.
- (2) Thermo-sensor-tests in stainless steel tubes may create false positive results. Small volumes of NCG heat up very quickly through the metal walls which have a 100-times higher heat transfer velocity in comparison to plastic or textiles.

3.37-55

154

gke

U. Kaiser

12/2012

Required steam quality for sterilizers

Standard	Maximum amount of NCG [ml] at 20°C at 1 bar allowed in 1 kg condensate	Refers to a maximum volume % in the steam volume at 121°C and 134°C
EN 285	35 ml (3,5%)	0,000027% (0,00027‰)

Calculation reference:

Gas	20°C 1 bar	100 °C 1 bar	121°C 2 bar	134°C 3 bar
1 kg steam [liter]	1244 l	1700 l	898 l	618 l
NCG [ml]	35 ml	48 ml	24 ml	17 ml

Potential risks during a fractionated vacuum steam sterilization process

1. Unsatisfactory air removal during the fractionated vacuum cycle
(remaining air in the sterilization chamber)
2. Leakages on door seals, valves and other devices
(air returns into the chamber after the last vacuum cycle)
3. Air transfer through the door seal, if the sealing is pneumatically
actuated
(if steam is used to pressurize, the sealing problem does not occur)
4. Non-condensable gases are introduced together with steam
(malfunction which is usually undetected during the sterilization cycle)
After sterilizer and eventually steam generator are switched off, non-
condensable gases develop in the pipes between steam generator and
sterilizer and in the steam generator, and get into the sterilizer after
starting again.

Non-condensable gases (NCG) in steam

Type of gas	Origin	Effect	Elimination
air ca. 80% N ₂ 20% O ₂	dissolved air in water (ca. 25 ml air per 1 l water)	The produced amount of air depends on the supply of injected water in the steam vessel. (air peaks after water injection)	Degasing of the injected water by heating up to 90 - 95°C before injecting into the steam vessel
air ca. 80% N ₂ 20% O ₂	In steam boilers and steam pipes before start-up	When not in operation steam boilers and steam pipes fill up with air	During start-up steam pipes need to be purged by a warm-up cycle to remove the air.
carbon dioxide CO ₂	water containing hydrogen-carbonates	During the heating process, CO ₂ and carbonate salt (white covering) are formed $\text{Ca}(\text{HCO}_3)_2 \rightarrow \text{CaCO}_3 \downarrow + \text{CO}_2 \uparrow + \text{H}_2\text{O}$ (Invisible for NCG detectors, because CO ₂ dissolves in condensate)	Degasing or de-mineralization by Ion-exchange or possibly together with Reversed Osmosis (RO) RO alone is <u>not</u> sufficient !
hydrogen H ₂	corrosion of metals	permanent small amounts of NCG and of flying rust in iron pipes (seldom)	Adjust pH-value > 7 in buffer solutions Remove chloride and other chelate complexing agents from feeding water
Superheated steam	pressure reduction in steam pipes	The superheated steam is unable to condensate until it is reaching its condensation point.	Install cooling line after pressure reduction.
	Hydratization of porous goods (heat generation by taking up water)		Don't dry porous goods with dry heat before the sterilization. Conditioning in normal humidity is required or moisten the goods before sterilization

3.41-55


158

U. Kaiser

07/2011

Recording of non-condensable gases (NCG) concentration in steam



3.42-55

159

U. Kaiser

07/2011

Water pre-treatment for steam generation (1)

1. Softening process using a cation exchanger

- Exchanges Mg^{2+} and Ca^{2+} against Na^+
- Regeneration by common salt (NaCl)

2. Reversed Osmosis

- retains salts but lets all gases, e.g. air, CO_2 , pass through the membrane

3. Mixed bed ion exchanger

- Exchanges all cations and anions against H^+ + OH^- (results in H_2O)

The water quality shall be measured with a conductivity meter, should not exceed 5 – 10 $\mu S/cm$ and should be monitored continuously in the CSSD. This treated feeding water still contains dissolved air.

4. Degassing before steam generator feeding

- heats up the water to 90 – 95°C
- degasses the water
- no loss of energy because heat up in the steam generator is anyway necessary later on

Different steam generators

1. Steam generation in the sterilizer

- directly in the chamber (small sterilizers)
- steam generator outside of the chamber, but short pipe to the chamber by electric heating

2. Decentral steam generation

- steam generator heating with oil, gas or electricity only for the sterilizer

3. Central steam generation

- central (expensive water preparation)


4. Secondary steam generation

- The primary steam is used as heating source
(The feeding water must be treated for the steam sterilizer)

Water pre-treatment for WDs and steam generation (1)

1. Water softening

in a reservoir with cation exchanger resin balls

Ca^{++} is exchanged against 2Na^+

3.47-55

331-1

gke


U. Kaiser

07/2012

Water pre-treatment for WDs and steam generation (2)

2. Reversed Osmosis

with plastic membrane (filter)

3.48-55

331-2

gke


U. Kaiser

07/2012

Water pre-treatment for WDs and steam generation (3)

3. Ion exchanger

in a reservoir with cation and anion exchanger resin balls, removes all inorganic salts from the feeding water for steam generation

3.49-55

332-1

gke


U. Kaiser

07/2012

Water pre-treatment for WDs and steam generation (4)

4. Degasser

removes the dissolved air

3.50-55

332-2


gke

U. Kaiser

07/2012

Measuring the conductivity

- Electrical conductivity is the ability of a material to conduct electricity
- Water conducts electricity, if salts (positive and negative ions) are dissolved in it.
- The more salt is dissolved, the higher is the conductivity. This is a measuring unit for the amount of dissolved salt in the water.
- The unit of conductivity is Siemens per meter (S/m) or mS/cm or $\mu\text{S/cm}$.
- $1 \text{ S/m} = 10 \text{ mS/cm} = 10.000 \mu\text{S/cm}$

Example for the conductivity of water:

- Absolute pure water : $0,05 \mu\text{S/cm}$ bis $0,1 \mu\text{S/cm}$
- Tap water: $300 \mu\text{S/cm}$ bis $1.000 \mu\text{S/cm} = 1 \text{ mS/cm}$
- Sea water: ca. $50.000 \mu\text{S/cm} = 50 \text{ mS/cm}$
- Demineralized water for
 - steam supply: $\leq 5 \mu\text{S/cm}$ (acc. EN 285)
 - end flushing: $< 15 \mu\text{S/cm}$ (acc. Guidelines of DGKH, DGSV, AKI)

3.51-55

587


gke

U. Kaiser

06/2012

Requirements for a successful sterilization process

1. Achievement of the kill kinetics, e.g.:
 $F_{121^\circ\text{C}} = 15 \text{ min overkill}$ in a steam sterilization process
2. Penetration of the sterilizing agent to all inner and outer surfaces

Test of the:

1. Kill kinetics by parametric release and/or biological and/or chemical indicators
2. Use of a process challenge device (PCD) to test the worst case penetration conditions using biological or chemical indicators or thermo-sensors inside the PCD.

3.52-55

162

gke

U. Kaiser

07/2005

**Pressure Dependence of Water Boiling Temperature;
Temperature Dependence of Heat Capacity and Steam Volume**

Boiling Temperature [°C]	Pressure [kpa] [bar]		Over-pressure [bar]	Heat Capacity [kcal/ kg steam]	Steam Volume [l/ kg steam]	Remarks
46	10	0,1	- 0,9		14560	Typical vacuum in a sterilizer
92	75	0,75	- 0,25		2200	Pressure on a mountain at a height of 2.500 m
100	100	1,0	0,0		1673	Standard pressure
121	200	2,0	1,0	647	885	Standard sterilization procedure at 121°C
134	300	3,0	2,0	659	606	Standard sterilization procedure at 134°C
145	400	4,0	3,0		446	

3.53-55

163

U. Kaiser

07/2005

**Reasons for wet goods
at the end of steam sterilization processes**

Reason	Remedy
Goods have been wet before packaging	Dry all goods to be sterilized before packaging
Condensate from the above goods is dropping down into the packages	Use completely covered metal container or cover each level of the sterilizer rack with a plate to let the condensate run outward
Wet steam, condensate aerosols enter the packages together with the steam	Use a pressure reducing valve before sterilizer and condensate traps to dry the steam
Condensate at the package/ container is separated from the goods and accumulates to a condensate "puddle" at the bottom of the container or soaks the bottom part of the packaging	Wrap the goods with a water-absorbing material (absorbent; cloth) to keep the condensate in the near of the goods and to use the heat of the goods to evaporate the condensate
At the end of the sterilization process steam remains in the package and condensates with gradual cooling	Use pressure-difference-process and flushing with air to remove steam at the end of the sterilization process

3.55-55

165


U. Kaiser

07/2005


Porosity of packaging materials

Woven packaging


Non-woven packaging

Paper

Paper stabilized with plastic fibers

- has large and small pores
 - when getting wet, paper loses its plain surface structure
 - cellulose fibers expand
 - after drying the paper does not get plain again
 - the pore size structure has changed
- plastic fibers stabilize the pore size distribution

13.1-2

309



gke

U. Kaiser

10/2010

Quality of packaging materials

Germs have a diameter of approx. $1\ \mu\text{m}$, pores $> 1\ \mu\text{m}$ may let germs pass

13.2-2

310

gke

U. Kaiser

10/2010

128.000 infections after surgeries each year

Quality survey in German CSSDs

Essen/ Münster, Germany. Alerting news in the Central Sterilization Service Department (CSSD) of hospitals: Due to insufficient sterilized surgical instruments the transfer of Creutzfeldt-Jakob disease (CJD) has been reported from England and Australia. For operations with a high risk of prions transfer, the usage of disposable instruments is mandatory.

In Germany on average experts assume the risk of 128000 post-operative infections. The affected patients stay for an additional 2 weeks in hospital causing additional costs in form of medications, medical personnel resources, loss of wages etc.

On invitation from the Centre for Hospital Management (CKM) of the Westfälische-Wilhelms-University in Münster, Germany 100 hospital managers, scientists, industrial government representatives discussed the risk of sterilization supply – in the area of conflict between legal requirements and cost efficiency.

The symposium director, Prof. Dr. Dr. Wilfried von Eiff, Director of CKM, outlined the different perspectives under which the CSSD in German hospitals has to be seen. „Zero fault quality“ is basically attributed to the controlled organizational process and the risk conscious behaviour of the hospital staff. Therefore the smallest faults have to be detected and eliminated.

From an economic point of view, the CSSD should be seen as an investment to avoid patient risks, rather than an expense factor. “In this area of conflict, an intelligent way between quality, risk and expenses has to be found”, mentioned von Eiff.

The government supervision of Nordrhein-Westfalen and Niedersachsen (2 Northern counties in Germany) analysed the sterilization equipment in hospitals and the quality of cleaning of endoscope units, e.g. at practising doctors. The preliminary survey results show that the requirements are not always met. More than 40% of hospitals which reprocess their instruments do not yet meet the legal requirements. 75% of these institutions can only reach the legal technical requirements by purchasing new sterilizers.

In many hospitals the CSSD is „the unpopular child in the basement“. Obviously the decision makers do not know what financial values exist in the CSSD: in a 400-bed hospital the value of instruments adds up to almost 5 Mio. €, added by 3 Mio € for service supply. Also the risks of incorrect work in the CSSD are not considered: incorrect trays and non-availability of instruments increase the costs by 200000 € per 10000 surgeries.

As a result of the CKM-Symposium a 16 point program has been developed, to support hospitals establishing a requirement oriented and efficient sterilization department.

Further information:


www.krankenhaus-management.de

Prof. Dr. Dr. Wilfried von Eiff,

University of Münster, Germany, Tel.: +49 (0) 2 51/83-3 14 40

Water Provision for Steam Generation in Steam Sterilisation

Practical tips from the Sterilisation Section of the German Society for Hospital Hygiene – DGKH e.V.

Extract from Annex 1 of the Recommendations by the German Society for Hospital Hygiene (DGKH) for validation of steam sterilisation processes revised by R. Fleischhack, H. Hahmann, P. Kober*, U. Kaiser, E. Drenthöfer

The quality of the water used as feed water is of special importance for the steam quality and hence for the sterilisation process. It is determined essentially by the method used to process the water. Water processing is intended as a means of fully or partially removing the non-condensable gases, hydrogen carbonates and any other dissolved mineral salts from the water. DIN EN 285 specifies corresponding limit values for the chemical and gaseous components permitted.

Information must be obtained from the regional water supply service on the composition of the drinking water (results of analyses) and on the processing methods employed.

How oxygen is supplied to remove iron and manganese may be of special interest in the context of steam generation for sterilisation. Compressed air and oxygen dosages are examples of the methods employed.

The former method, especially, transports a large quantity of dissolved gases (in particular oxygen and nitrogen) to the consumer. Special attention must be paid to this water when subjecting it to further processing to obtain feed water for steam generation.

The customary method of processing water for Central Medical Devices Processing (CMDP, CSSD) in the hospital is as follows:

- Water softening
- Reverse osmosis
- Ion exchangers/mixed-bed exchangers
- Storage containers for demineralised water
- Conveyance to the consumers in CMDP, CSSD

Water softening removes the hydrogen carbonates present in the water as they are broken down by heat to CO_2 and carbonates (limescale).

The reverse osmosis system is able to remove around 90% of the salts from raw water but not the CO_2 , oxygen and nitrogen dissolved in the water.

Postconnected ion exchangers (series connection of 2 recommended) can remove any remaining salts as well as the CO_2 . Silicates, however, pass through ion exchangers before onset of a change in conductance. This is counteracted by, as proposed, connecting the ion exchangers in series. Large quantities of dissolved nitrogen are not trapped by any of the aforementioned processing aids. Hence special measures must be taken here.

Impeccable functioning of the steriliser can be assumed for conductance values that are continually below 5 mS/cm in the feed water used for the steam generator. In line with recent findings, this value deviates from that specified hitherto in DIN EN 285, Annex B.

Processed water has a very high affinity for all gases present in its immediate environment and thus quickly dissolves them back. This occurs especially if the water surface in the storage container is greatly enlarged due to the collision resulting from the addition of water. This tendency to dissolve back gases must therefore be counteracted by observing what is known as "subsurface filling". This ensures that the surface of the water in the storage container is kept relatively calm and that not too much air, or especially carbon dioxide, is dissolved back.

When the water, which has been processed and stored in a storage container, is supplied to the consumer it contains dissolved oxygen and nitrogen (approx. 25 – 30 ml gas/l water). This proportion has no adverse effect on using demineralised water in washer-disinfectors.

To generate the sterilant "saturated steam", the presence of such gaseous

additions in the form of non-condensable gases (NCGs) can have important implications. Hence the latter should be removed by as far as possible resorting to thermal degassing. Otherwise this quantity of dissolved gas would be conveyed into the chamber together with the steam on each subsequent occasion that water is supplied.

The conductance value of the first ion exchanger should be measured and documented daily as a proof of the conductance of the feed water mediated by carbon dioxide. A limit-value signalling device should be used to round off conductance measurement.

Conductance mediated by carbon dioxide will be eliminated if processing is properly effected using postconnected ion exchangers and thermal degassing of the feed water at a boiling point temperature for water under standard conditions. In this respect, the conductance value measured is an indicator of the steam quality.

Ion exchangers have a finite capacity and must be replenished from time to time. The silicate problems mentioned can be overcome by connecting the ion exchangers in series as recommended.

If there are noticeable signs of adverse effects on the sterile supplies and chamber (e.g. discolorations, deposits) following sterilisation, carrying out water analysis as per DIN EN 285 is advisable, in particular before the initial validation, regardless of the type of steam supply. If this water analysis is unable to identify this problem, other troubleshooting measures must be taken (e.g. contamination of the feed water storage container, use of inhibitors, rustproof pipes, etc.).

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Sterilization of minimal invasive surgical (MIS-) instruments with splits, screw threads and sealing gaskets using different maintenance, protection and lubricating agents in steam sterilization processes

Necessary conditions for a successful steam sterilization processes

- All surfaces to be sterilized have to be wet.
- Temperature-time-integral of water minimum 121°C, 15 min ($F_0 \geq 15$ min) at all surfaces to be sterilized is required.

The sterilizing agent is water which has achieved a temperature-time-integral.

The condensation of steam is no requirement. However, it allows a very quick heat transfer.

There is no sterilization in steam sterilization processes possible, if:

- steam cannot condensate to water due to overheating (e.g. hygroscopic condensation in cellulose fibres, overheated steam etc.)
- the surfaces to be sterilized are sealed and thus the water cannot get to the surfaces of the instruments, e.g. because of:
 - bio-films
 - maintenance and protection agents
 - lubricants
 - sealing gaskets
 - accumulation of non-condensable gases (NCG) in partially closed spaces


Complex instruments:

- have tight fitted surfaces
- have screw threads
- have elastic sealing gaskets
- have sliding sealing gaskets (e.g. valves at trocars)
- are surface covered with maintenance and protection agents and/or lubricants
- have narrow lumina [*]

These versions have been tested in a research project.

[*] see publication : Kaiser U, Gömann J: Investigation of Air Removal from Hollow Devices in Steam Sterilization Processes. ZentrSteril 1998; 6 (6): 401-413

Two stainless steel plates are contaminated with *G. stearothersophilus* suspension, wrapped and sterilized in a steam sterilization process (134°C, 5 min)

4.6-13

171

gke

U. Kaiser

01/2006

Test in a split of inoculated stainless steel plates

In the middle the plates have been inoculated in a circle of 5 cm diameter with:

100 µl suspension
G. stearothersophilus

Population:
1.4 x 10⁸ CFU/ ml
D₁₂₁-value = 2.1 min
z-value = 10°C
F_{BIO}-value = 15 min

4.7-13


172

gke

U. Kaiser

07/2005

**Test with pressed stainless steel plates
with an inserted elastic silicone sealing gasket**

4.8-13

173

gke

U. Kaiser

07/2005

Test with a M8 screw thread of 20 mm length

4.9-13

174

gke

U. Kaiser

07/2005

Different maintenance and protection agents and lubricants

- Validation according EN ISO 17664 required by the manufacturer.
- The manufacturer has to provide directions for use and a validation test report according EN ISO 17664 on request.

4.11-13


176


U. Kaiser

09/2009

Incubation of a stainless steel plate in culture media

+ growth = NON-STERILE

- No growth = STERILE

4.12-13

177

U. Kaiser

09/2009

Results of the microbiological tests in steam sterilization processes

Sterilization procedure	Gravity cycle		Process with fractionated vacuum [**]	
	Stainless steel plates screwed up	Screw thread	Stainless steel plates screwed up	Screw thread
- sterilized	+	-	-	-
- sealed with silicone plates and sterilized	+	n.a.	+	n.a.
- treated with maintenance and protection spray based on aqueous suspension and sterilized	+	-	-	-
- treated with maintenance and protection oil [*] and sterilized	+	-	-	-
- treated with lubricating Grease and sterilized	+	+	+	+

- No growth = **STERILE** + Growth = **UNSTERILE**

* Those maintenance and protection oils containing detergents can absorb water. Lubricating grease without detergents does not contain water.

** B1-cycle according to EN ISO 11140-4

Results


1. Narrow splits which are not sealed, can be sterilized in fractionated vacuum cycles but not in gravity cycles.
2. Sealed surfaces cannot be sterilized in any steam sterilization process.
3. Maintenance and protection sprays which consist of aqueous suspensions or non polar protection oils, which contain detergents and can absorb water, do not hinder sterilization.
4. Protection oils and lubricating grease which cannot absorb water and sealed surfaces inhibit sterilization.

**Air removal from hollow goods
(for example: catheters, MIS-Instruments etc.)**

**Tube system both sides open,
for example 3 m length**

Indicator strips to monitor steam penetration

**Tube system half length one side closed,
for example 1,5 m length**

In steam-sterilization processes both systems have the same difficulty of air removal.

5.1-13

179

gke

U. Kaiser

07/2005

**Examination of air removal and steam penetration in tubes
depending on the tube dimensions and the air removal process**

Tube diameter [mm]	Tube length [m]	Colour change of the chemical indicator at the following number of vacuum cycles from 100 to 950 mB					
		1	2	3	4	5	10
1	0,5	0	100	100	100	100	100
	1	0	100	100	100	100	100
	2	0	90	100	100	100	100
	3	0	90	100	100	100	100
	4,5	0	75	100	100	100	100
2	0,5	0	75	100	100	100	100
	1	0	75	100	100	100	100
	2	0	30	100	100	100	100
	3	0	10	80	100	100	100
	4,5	0	20	90	100	100	100
3	0,5	0	80	100	100	100	100
	1	0	30	100	100	100	100
	2	0	25	100	100	100	100
	3	0	0	25	100	100	100
	4,5	0	0	50	90	100	100
4	0,5	0	50	100	100	100	100
	1	0	0	100	100	100	100
	2	0	0	100	100	100	100
	3	0	0	90	100	100	100
	4,5	0	0	20	60	90	90
5	0,5	0	20	100	100	100	100
	1	0	0	100	100	100	100
	2	0	0	100	100	100	100
	3	0	0	10	90	100	100
	4,5	0	0	0	10	10	25

Construction of the process-challenge device (PCD)

Test device: **gke**-brass test device

(batch monitoring system)

Tube: PTFE-tube, various lengths and diameters

Air removal

minima and maxima

of vacuum cycles: 100 to 950 mB


Number of

vacuum cycles: various (1 to 10)

Rate of pressure rise: 1000 mB ± 200 mB/ min
(Sterilization at 134°C for 3:30 min)

Results

Increasing length and diameter of the tube causes more difficult air removal of the PCD.

5.5-13

183

gke

U. Kaiser

05/2012

Examination of air removal and steam penetration in tubes depending on the tube dimensions and the air removal process

HPR L x d [cm ²]	Tube diameter [mm]	Tube length [m]	Colour change of the chemical indicator at the following number of vacuum cycles from 100 to 950 mB						
			1	2	3	4	5	10	
5	1	0,5	0	90	100	100	100	100	
10	1	1	0	90	100	100	100	100	
10	2	0,5	0	75	100	100	100	100	
15	3	0,5	0	50	100	100	100	100	
20	1	2	0	50	100	100	100	100	
20	2	1	0	50	100	100	100	100	
20	4	0,5	0	40	100	100	100	100	
25	5	0,5	0	20	100	100	100	100	
30	1	3	0	30	100	100	100	100	
30	3	1	0	25	90	100	100	100	
40	2	2	0	25	90	100	100	100	
40	4	1	0	10	50	100	100	100	
45	1	4,5	0	20	90	100	100	100	
50	5	1	0	0	10	100	100	100	
60	2	3	0	0	30	75	100	100	
60	3	2	0	0	30	50	100	100	
80	4	2	0	0	20	60	100	100	
90	2	4,5	0	0	10	50	100	100	
90	3	3	0	0	0	40	40	100	
100	5	2	0	0	0	10	25	100	
120	4	3	0	0	0	20	75	100	
135	3	4,5	0	0	0	0	25	100	
150	5	3	0	0	0	0	10	100	
180	4	4,5	0	0	0	0	0	20	
225	5	4,5	0	0	0	0	0	10	

Table of measured values according to the increasing HPR

**Construction
of the process-challenge-device (PCD)**
Test device: **gke** brass test device
 (batch monitoring system)
Tube: PTFE-tube, various lengths
 and diameters

**Air removal
minima and maxima
of vacuum cycles:** 100 to 950 mB

**Number of
vacuum cycles:** various (1 to 10)

Rate of pressure rise: 1000 mB ± 200 mB/ min
 (Sterilization at 134°C for 3:30 min)

Results
 The same HPR-value of different PCD's shows
 approx. the same result, for example:

$$1 \times 4 = 40$$

$$4 \times 1 = 40$$

$$2 \times 2 = 40$$

5.6-13

184

U. Kaiser

07/2005

Dependence of the inner tube surface and the tube volume on the diameter of the tube

Parameter [Dimension]		
Diameter [mm]	2	4
Length [m]	1	1
Inner volume [ml]	3,1	12,6
Inner surface [cm ²]	62	126
Ratio Volume : Surface [ml/ cm ² = cm]	1 : 20	1 : 10

5.12-13

190

U. Kaiser

07/2005

**Requirements for sterilization monitoring
in steam sterilizers with fractionated vacuum**

1. Validation / revalidation (annually)

Test of the sterilizer at the manufacturer, after installation and substantial repairs = operation and installation qualification (OQ & IQ) = Commissioning
Performance qualification of all processes / programs using worst-case load configurations = performance qualification (PQ)

2. Function test of the sterilizer (during each start-up)

Function test at start-up of the sterilizer to prove air removal and steam penetration using a BD- and vacuum test.

3. Routine monitoring of each batch (in each cycle)

Recording the basic technical parameters (for steam sterilizers: pressure and temperature versus time and monitoring of successful steam penetration)

4. Documentation using a quality management system (QMS)

Sterilized goods must be released by trained personnel according to defined standard operation procedures (SOP).

**Robert Koch Institute (RKI) Recommendation
Requirements to the Hygiene during
Reprocessing of Medical Devices, 10/2012**

Quote of chapter 1.4

„Securing the quality of the reprocessing processes to be used“
(page 1251)

The quality of the reprocessing is secured in dependence of the individual cleaning, disinfection and sterilization procedures by

- a) validation (consisting of installation [IQ], operation [OQ] and performance qualification [PQ]),*
- b) regular routine monitoring (e.g. daily),*
- c) batch-related routine monitoring,*
- d) monitoring and checking of the process parameters,*
- e) maintenance, calibration, if necessary adjusting, servicing and*
- f) regular process validations (revalidation)*
- g) event-related validation (performance qualification because of special reasons)*

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1. Validation

EN ISO 17665-1 describes the validation and routine monitoring of steam sterilization processes. It is valid for all processes used in industry, laboratories and health care institutions, up to disinfection of liquids and waste. Therefore the requirements of the standard are very general, since depending on the sterilization process (e.g. with or without fractionated vacuum) and the load configuration (e.g. instruments, liquids, waste) different monitoring methods have to be used.

Therefore the general description of the validation and monitoring procedure in the standard always has to be adapted to process and load configuration.

2. Function test of the sterilizer

Function test of the sterilizer at start-up

Clause 12 of the validation standard EN ISO 17665-1
(Maintenance process effectiveness) quotes:

12.1.6: „If the sterilization process relies on the removal of air from the sterilizer chamber in order to achieve rapid and even penetration of steam into the sterilizer load, a steam penetration test shall be carried out each day before the sterilizer is used.”

The standard does not prescribe a special test for routine monitoring, therefore usually two type tests described in the sterilizer standard EN 285 are used:

- Bowie-Dick-Test according to EN 285 Part 17
- Hollow Load Test according to EN 867-5 and EN 285

gke offers a test device for the BD- and the hollow load test simultaneously.

3. Routine monitoring of each batch

Routine monitoring of each batch

In the validation standard EN ISO 17665-1 the following text is quoted:

10.1. „Routine monitoring and control shall be performed on each operating cycle.“ (each cycle = each batch)

10.3. „Steam penetration“

Clause 11: Product release from sterilization

11.1. „Procedures for the [...] product release from the sterilization process shall be specified.“

This specification must define the condition that the sterilization process is conforming to specification (see 9.5.2, **10.3**, 10.5 and 10.6). If one condition is not fulfilled, the product must be marked as faulty.

3. Routine monitoring of each batch

Therefore the standard EN ISO 17665-1 requests for routine monitoring:

- All critical parameters of each batch must be monitored and tested for sterility.
- The parameter „steam penetration“ according 10.3 is depending on the complexity of the load, one of the most critical parameters and must be tested in each batch.

gke offers for this monitoring special load-related PCD systems. These PCD systems must be adapted to the requirements of steam penetration in the load.

- Monitoring of steam penetration with air detectors is in instruments neither qualitatively nor quantitatively possible. The amount of non-condensable gases in steam is only a part of the steam penetration problem, the correct air removal of individually different complex instruments cannot be measured with air detectors.

4. Documentation using a quality management system (QMS)

Information and documents, which must be available for release:


Batch printout of the sterilizer with pressure and temperature over time

Process indicators – provide only information that the package has passed the sterilization process. No information about sterility.

Prove of sufficient steam penetration (depending on the load)

- a. Use of biological or chemical indicators in packs to prove, if solid instruments are sterilized
- b. Use of Indicator systems in PCDs (Medical Device Simulator = MDS or Batch Monitoring System = BMS), if hollow instruments, complex loads or waste are sterilized
- c. Use of biological indicator suspension and/or temperature sensors if liquids are sterilized

Compliance with release conditions must be batch-related documented.





**Different types of the Bowie-Dick-Test
Comparison Europe - USA**

Country	Standard	Size [cm]	Weight [kg]	BD-Test-program inside the sterilizer	Required method of simulation test		
					Standard for test	name	Test method
Europe	EN 285, part 17	25 x 35 x 20	7 kg ± 10%	134°C, 3.5 min or 121°C, 15 min	EN 867-4 = EN-ISO 11140-4	Steam penetration test	- Air removal - Leaks - Non-condensable gases
USA	AAMI/ANSI ST 79	ca. 24 x 35 x 29	4 kg ± 200g (± 5%)	132°C, 3 min	ISO 11140-5	Air removal test	- Air removal

The American test package has approximately 1/2 of the weight of the European test package and is therefore less sensitive in testing the air removal and steam penetration.


6.10-32 194 *gke* U. Kaiser 06/2010

Air removal of porous and hollow goods using different steam sterilization processes			
BD-Test package according to EN 285	Hollow load test according to EN 867-5	+ - = penetration successful = penetration not successful	
+	+		Sub-atmospheric cycle according to EN 867-4, B.1  Pressure [mbar] Time [min]
pass	pass		
+	-		Trans-atmospheric cycle according to EN 867-4, B.3 with a sterilization temperature of 134°C  Pressure [mbar] Time [min]
pass	fail		
+	-		Trans-atmospheric cycle according to EN 867-4, B.3 with a sterilization temperature of 121°C  Pressure [mbar] Time [min]
pass	fail		
From the publication Gömann/ Menzel/ Kaiser – Central Sterilization 3/2001			
During the modification of the European Standard EN 285 for big sterilizers the hollow load PCD according to EN 867-5 was included as standard-test additionally to the BD-Test pack.			
6.11-32	195		U. Kaiser 07/2009


Technical principle of the *gke* Compact-PCD™

This part is integrated into the inner part of the plastic case

Volume		
Detector	Stainless steel tube 1m	Plastic housing
0.2 ml	3 ml	100 ml

6.13-32

197

gke

U. Kaiser

10/2010

Why Batch Monitoring?

Why using Batch Monitoring Systems? (1)

- After validation batch monitoring is not necessary, if the process remains stable.
- A validated process must be 100% reproducible over a longer period.
- Not all critical parameters of the cleaning and sterilization process are 100% reproducible (different soils, load, different steam quality, potential leaks, etc.).
- Therefore 100% reproducibility is not guaranteed.
- This process risk has to be compensated checking all the critical non-reproducible parameters which are not controlled by the equipment with a suitable monitoring system.

6.16-32

358

gke

U. Kaiser

05/2012

Why Batch Monitoring?

Why using Batch Monitoring Systems? (2)


Batch monitoring systems are used, if a device, i.e. sterilizer or washer/disinfector, does not monitor and document all critical variables.

A batch monitoring system can monitor one, several or all critical variables, so that the result of a batch monitoring system allows the release of a load (i.e. in a steam sterilization process monitoring of air removal, temperature-time-integral (F_0 -value) and the condensation of steam to water) is sufficient.



It is a necessary requirement that the batch monitoring system is adjusted according to the load. Therefore the batch monitoring system has to be validated with suitable methods, i.e. using the method described in DIN 58921 which is a new work item in the ISO/TC 198 technical committee to become an ISO standard.

If the load consists of different products, the validation of the batch monitoring system must relate to the most difficult load configuration for a cleaning or sterilization process.

Limitations of biological and chemical indicators at different positions in steam sterilization cycles

Process Challenge Devices (PCDs) for different applications

Reference	is simulated by	Process Challenge Device
Medical Device 	 Using the simulation test procedure of the standard DIN 58921	Medical Device Simulator 
Batch  Defined Load Configuration		Batch Monitoring System 

6.19-32

200

J. Metzing

12/2010

Class 2 indicators for type tests, MDS¹ and BMS²

MDS or BMS test-system are used to monitor sterility of a medical instrument or load in sterilization processes:

Process
Challenge + Indicator system = **I**ndicator
Device

An indicator is an object in its final form in which it is intended to be used.
(Definition in EN-ISO 11140-1 for class 2 indicators)

¹ MDS = Medical Device Simulator

² BMS = Batch Monitoring System

6.20-32

201

U. Kaiser


03/2007

**Robert Koch Institute (RKI) Recommendation
Requirements to the Hygiene during
Reprocessing of Medical Devices, 10/2012**


Quote in attachment to annex 4
„Special requirements for the use of chemical indicators“
(page 1274)

PCDs have been developed to simulate a special difficulty of steam penetration in a sterilization process [...] There is no universal PCD to be used for all processes and applications.

**Test method to check if a Medical Device Simulator (MDS)
or Batch Monitoring System (BMS)
is equivalent to a Medical Device (MD) or defined load
according DIN 58921**

Conceptual approaches for routine monitoring with a batch monitoring system (BMS)

6.24-32


204

gke

U. Kaiser

07/2009

Use of a batch monitoring system (BMS)

6.25-32

205

gke

U. Kaiser

02/2011

**Different types of monitoring systems
testing the efficacy of steam sterilization processes**

	Type of test	Monitoring potential	Limitations
1	Parametric release: recording the temperature- pressure-time window	Temperature-time- windows checking the kill of germs	No ability to test the chamber for non-condensable gases in steam coming from various sources
2	Biological indicators		Test only at the place where indicators are located
3	Non-biological all parameter integrating indicator (emulates biological indicators according to EN 867-1 class D)		
4	Bowie-Dick-Test	Steam penetration test for porous loads and solid instruments	No steam penetration test for complex hollow devices
5	Specific process challenge device (PCD) simulating complex hollow devices according to EN 867-5	Steam penetration test of solid instruments, porous loads and hollow devices and kill kinetics	Comparison is necessary between the penetration characteristics of the hollow devices to be tested and the PCD used

6.26-32

206

U. Kaiser

07/2005

**Robert Koch Institute (RKI) Recommendation
Requirements to the Hygiene during
Reprocessing of Medical Devices, 10/2012**

Quote in attachment of annex 4
„Information for the use of biological and chemical indicators“
(page 1273)

*A successful validation does not secure that the sterilization conditions
are always kept during routine operation. [...]*

*Therefore suitable routine monitoring must be carried out
(see EN ISO 14937, clause 10), to secure that the validated process is
always carried out.*

6.27-32

364

J. Metzger

02/2013

**Robert Koch Institute (RKI) Recommendation
Requirements to the Hygiene during
Reprocessing of Medical Devices, 10/2012**

Quote of checklist
„Implementation and operation of small sterilizers
to reprocess medical devices“
Heading „Batch-related tests“
(page 1272)

Test of the complete and correct process run

- *Test and documentation of the result of the process indicators (class 1)*
- *Test and documentation of the process parameters*
- *Test and documentation of the result of the indicator*
 - *Critical A instruments: without PCD (class 5)*
 - *Critical B instruments: with PCD, e.g. helix test (class 2)*

**Robert Koch Institute (RKI) Recommendation
Requirements to the Hygiene during
Reprocessing of Medical Devices, 10/2012**

Quote in attachment of annex 4
„Special requirements for the use of chemical indicators“
(page 1274)

The level of difficulty of class 2 indicators is depending on the combination of the chemical indicator and the PCD components. Each variation, e.g. the use of another indicator, can cause that malfunctions are not detected. Therefore only components defined from the manufacturer can be combined. This applies also for systems with a reusable PCD.

Principles of non condensable Gas (NCG) Detectors

NCG-detectors are offered with the argument, that the necessary test for steam penetration of each batch according to EN ISO 17665-1 clause 11.1 is fulfilled by these detectors. This information is:

- technically incorrect and
- not in line with the content of the standards.

The use of NCG-detectors requires the knowledge of the different technical solutions and requirements of the standards.

8.1-15

315


gke
GAS DETECTION SYSTEMS

J. Metzging

08/2010

Principles of non condensable Gas (NCG) Detectors

(1) Integral Method according EN 285

- (1) steam is split and condensed in a heat exchanger
- (2) water is acidified with H_3PO_4 , that no CO_2 dissolves and can be detected
- (3) Condensate and NCG are collected in a transparent cylinder
- (4) The condensate amount is measured with a balance

Advantages

- NCG are collected over time and the volume integrates
- if acid is added, CO_2 is detectable

Disadvantages

- the method integrates NCG over time, therefore no NCG maxima are detectable
- correct results are only possible, if a homogenous steam/NCG current is in the steam pipe
- no information available when NCG are entering the process, NCG during the air removal and/or plateau phase are not critical for the process (but may show fail results in a pass cycle)

¹⁾ EN 285 requires from a sterilizer still to operate with a maximum concentration of 35 ml NCG/1 kg condensate.

8.2-15

255


gke
GAS DETECTION SYSTEMS

U. Kaiser

08/2010

Principles of non condensable Gas (NCG) Detectors

(2) Differential Method

- (1) steam is split and condensed in a heat exchanger
- (2) condensate and possibly NCG pass a cuvette with a light beam
- (3) light beam is deflected differently in condensate and NCG
- (4) Condensate amount is online measured

Advantages	Disadvantages
- NCG-peaks can be detected over the time	- gas detector (GD) may only be validated above a precision of 5%. Quite often the GD is only validated against a PCD system, without information about the precision - CO₂-gas is undetectable because it gets dissolved again in the condensate - correct results are only possible, if a homogenous steam/NCG current is in the steam pipe. - no information available when NCG are entering the process, NCG during the air removal and/or plateau phase are not critical for the process (but may show fail results in a pass cycle)

8.3-15

256


gke

U. Kaiser

08/2010

Principles of non condensable Gas (NCG) Detectors

(3) Integral Method similar to EN 285

- (1) steam is split and condensed in a heat exchanger
- (2) condensate and NCG are collected in a transparent cylinder
- (3) light beam is deflected differently in condensate and NCG and the amount of NCG is detectable in the cylinder
- (4) Condensate amount is online measured

Advantages	Disadvantages
- over the time NCG-peaks can be detected	- gas detector (GD) may only be validated above a precision of 5%. Quite often the GD is only validated against a PCD system, without information about the precision - CO₂-gas is undetectable because it gets dissolved in the condensate - correct results are only possible, if a homogenous steam/NCG current is in the steam pipe - no information available at what time NCG are entering the process, NCG during the air removal and/or plateau phase are not critical for the process (but may show fail results in a pass cycle)

8.4-15

257


gke

U. Kaiser

08/2010

Principles of non condensable Gas (NCG) Detectors

(4) Integral Method with pressure difference method

- (1) steam is split and condensed in a heat exchanger
- (2) Condensate is collected in tube of 75 cm length
- (3) Pressure is measured in the filled tube containing a condensate / NCG mix
- (4) NCG is measured by pressure depression compared to a tube filled with water

Advantages

- over the time NCG-peaks can be detected

Disadvantages

- gas detector (GD) may only be validated above a precision of 5%. Quite often the GD is only validated against a PCD system, without information about the precision
- CO₂-gas is undetectable because it gets dissolved in the condensate
- correct results are only possible, if a homogenous steam/NCG current is in the steam pipe
- no information available at what time NCG are entering the process, NCG during the air removal and/or plateau phase are not critical for the process (but may show fail results in a pass cycle)

8.5-15

328


gke

U. Kaiser

08/2010

Principles of non condensable Gas (NCG) Detectors

(5) Tests by Temperature Gradients

- (1) non-condensable gases may accumulate at cold places
- (2) between the non-condensable-gas-nest and the chamber a temperature difference may built up which can be recorded over time
- (3) the temperature difference allows a qualitative evaluation of non-condensable gases

Advantages

- simple electrical recording of the temperature difference vs. time
- not expensive

Disadvantages

- no quantitative method to test non-condensable gases (There are only indirect assumptions about the amount of non-condensable gases possible.)
- regarding the security of a sterilization process the interpretation of the results is not possible without additional parameters.
- no validation of the non-condensable-gas-test is possible.
- no relation to the load to be sterilized
- no information available when NCG are entering the process, NCG during the air removal and/or plateau phase are not critical for the process (but may show fail results in a pass cycle)

8.6-15

258

gke

U. Kaiser

08/2010

Principles of non condensable Gas (NCG) Detectors

NCG-Detectors:

- cannot be validated.
- do not detect CO₂ as NCG, because this gas dissolves in water.
- monitor the steam quality in one location but not the steam penetration into hollow instruments.
- take steam samples from the steam pipes or from the sterilizer chamber but not from inside the package in which NCG are accumulated by steam consumption and condensation.
- are not calibrated in their sensitivity to the requirements of the each load.
- do not only measure during the critical heat-up time but during the whole program and can show fail results of a pass cycle.

Principles of non condensable Gas (NCG) Detectors

1. NCG-Detectors cannot be validated

Test systems must be validated against a reference. There is no known method for the validation of an NCG-Detector for the NCG-concentration.

The sensitivity of NCG detectors does not record absolute concentrations of NCG but is calibrated using other test methods. If the detector signals a fail (depending on the fail-condition, see point 6), such a test is from a technological point of view incorrect.

Many NCG detectors use an integral measurement – a summary of NCG over a time interval.

Only few NCG detectors use a differential measurement. The value recorded is not synchronized with the sterilization parameters.

Principles of non condensable Gas (NCG) Detectors

2. NCG-Detectors do not detect CO₂ as an NCG, because dissolving in water

All detectors based on the assumption that steam condenses but NCG do not condense, can only measure those NCG not dissolving in water.

Gases like CO₂ dissolving in water (condensate) are not detected.

Therefore NCG-Detectors do not detect CO₂ and small amounts of O₂.

Risk of a false-positive test result!

Principles of non condensable Gas (NCG) Detectors

3. NCG-Detectors monitor the steam quality but not the steam penetration

The validation standard EN ISO 17665-1 requires successful monitoring of steam penetration into each load also into hollow instruments. Successful steam penetration is assured using a batch monitoring system consisting of PCD and indicator system.

An NCG detector measures the steam quality, e.g. coming from a steam pipe or sterilizer chamber. The conclusion that acceptable steam penetration into a hollow load is automatically assured, is wrong.

If air removal is insufficient, air remains in the load (i.e. MIC-instruments). Even if the steam quality is acceptable, steam penetration into hollow devices is blocked.

Risk of false-positive test results!

Principles of non condensable Gas (NCG) Detectors

4. NCG-Detectors take out steam samples from the steam pipes or from the sterilizer chamber but not from the package

Monitoring of steam penetration must be done inside of the load or rather inside of hollow instruments.

If NCG come into a chamber through leaks and the NCG detector takes samples out of a steam pipe, the failure cannot be detected.

Therefore many reasons for unsuccessful sterilization are not detected by an NCG-detector.

Risk of false-positive test results!

Principles of non condensable Gas (NCG) Detectors

5. NCG-Detectors are not adapted in their sensitivity to the requirements of a load configuration

A BMS has to be adapted to the specifications of a worst-case load configuration. This requirement is guaranteed by the selection of a suitable PCD validated according DIN 58921.

An NCG-detector cannot be adapted to different loads. It is a device that monitors the function of a sterilizer but does not provide any release criterion for successful sterilization of a load.

Principles of non condensable Gas (NCG) Detectors

6. NCG-Detectors do not only measure during the critical heat-up time but during the whole program


The most frequent reason for failures in steam sterilization processes is the introduction of NCG with steam. The NCG concentration in steam is not constant but may change heavily during the day.

The introduction of NCG is only critical during the heat-up time. During other times, e.g. during the fractionating or plateau phase, NCG do not jeopardize the process, because during this time there is no steam consumption in the load.

Risk of false-negative test results!

If these connections are unknown, there is the risk that NCG-detectors are considered higher than other test systems.

Principles of non condensable Gas (NCG) Detectors

To demonstrate the superior function of an NCG-detector NCG (air) is injected in the sterilizer during an uncritical time. This test should demonstrate the high sensitivity of the NCG-detector.

The detector, measuring steam in the chamber during the whole time, indicates this air injection. However the success of the sterilization is not jeopardized by this injection because the largest part of steam consumption in the load has been already consumed. MDS and BMS do not show any failure because the heat-up phase of the good is already nearly reached and NCG don't have a negative influence anymore.

During that comparison the NCG-detector will show a fail, the MDS or BMS will show a pass.

Differences between the usage of biological and chemical indicators (1/4)

In USA biological indicators (BIs) are the „Golden“ Standard for sterilization monitoring. The reason is that alive germs (BIs) similar like pathogenic germs are used to monitor the efficacy of sterilization procedures.

In most of the sterilization processes (except H₂O₂) it is very well-known which physical parameters are relevant to inactivate pathogenic germs and biological indicators.

If the parameters inactivating alive bacteria (i.e. in steam: temperature, time, water) are absolutely known, a biological indicator to prove kill is not anymore necessary, since those critical parameters can be controlled by other methods more accurately.

In steam sterilization temperature and time can be monitored by recording those physical parameters by the sterilizer, the existence of water at all surface to achieve sterility, must be checked, however at the most difficult penetration location in the load, independent if a biological indicator or any other test method for water is used.

In steam sterilization processes today chemical indicators are available which only react to the pass conditions, if the critical parameters temperature, time and water are achieved.

Therefore from a risk management point of view there is conceptually no difference using biological or chemical indicators.

14.1-6

335

U. Kaiser

01/2011

Differences between the usage of biological and chemical indicators (2/4)

Advantages of biological indicators (BIs):

- Direct inoculation inside of complex instruments possible
- Required in many low temperature sterilization processes, since adequate chemical indicators are not available

Disadvantages of BIs in steam sterilization processes:

- The indicator has to be removed from the package to be incubated, there is no data available after opening the pack
- Incubation time required, no immediate information about the process available, delay in time until information is available
- Rapid development of BIs has the deficiency in probability:
longer incubation time is increasing the probability to get correct results but also longer time is necessary
- In processes at 132 – 134°C BIs do not provide any data how long the process was running, since BIs are inactivated in less than 1 minute
- Cost including inoculation is 3 – 5 times more expensive than chemical indicators

14.2-6

336

U. Kaiser

01/2011

Differences between the usage of biological and chemical indicators (3/4)

Advantages of chemical indicators (CIs) as long as they conform with the requirements according EN ISO 11140-1 testing all critical parameters (class 5 + 6):

- Equivalent or better in pass/fail results in comparison to BIs
- Immediate response to all critical parameters
- Data is available at the end of the process
- Cost-effective in comparison to BIs
- Indicator can be stored for documentation
- No incubation necessary
- Stated values can be adapted according to the sterilization process used, which is not possible for biological indicators at 132 - 134°C

Disadvantages of CIs:

- No direct inoculation at worst-case locations possible
(but only needed for validation, not for routine monitoring)
- A lot of CIs are on the market without required specifications discrediting CIs,
(i.e. several CIs for steam sterilization pass in a dry heat process as well which should not happen)

14.3-6

337

U. Kaiser

01/2011

Pass/fail conditions of biological and chemical indicators


Example of pass/fail conditions of steam sterilization indicators with a stated value (SV)
at 121°C

Indicator Type	Standard	Fail-Time [min]	Pass-Time [min]
Biological indicators	EN ISO 11138-1 + 3	4,5	13,5
Chemical indicators class 5	EN ISO 11140-1	10,5	16,5
Chemical indicators class 6		15,5	16,5

Biological indicators

Chemical indicators class 5

Chemical indicators class 6

- fail
- pass
- uncertain (pass or fail)

14.4-6

115

U. Kaiser

11/2010

Pass/fail conditions of biological and chemical indicators

Example of pass/fail conditions of steam sterilization indicators with a stated value (SV) at 134°C

Indicator Type	Standard	Fail-Time [min]	Pass-Time [min]
Biological indicators	EN ISO 11138-1 + 3	0,33	1
Chemical indicators class 5	EN ISO 11140-1	2,25	3
Chemical indicators class 6		2,7	3

Biological indicators

Chemical indicators class 5

Chemical indicators class 6

■ fail
■ pass
■ ■ uncertain (pass or fail)

14.5-6

339

gke

U. Kaiser

01/2011

Differences between the usage of biological and chemical indicators (4/4)

History and current situation of biological indicators (BIs) in different countries

BIs were the first test method to prove sterility. Most hygienists and microbiologists are not updated and therefore still believe that only the biological indicator is a proof for sterility (reason for the "Golden" Standard).

In most countries world-wide including USA biological indicators are not only recommended but are mandatory to be used due to missing local education.

In the West European countries (Scandinavia, United Kingdom, Germany, Austria and Switzerland) biological indicators are only rarely used for routine monitoring in steam sterilization processes, only used for validation and in low temperature sterilization processes.

Update in know-how is necessary to opinion leaders world-wide, the current problem is, that in universities world-wide there is no research and education, therefore historical status is remaining.

14.6-6

338

gke

U. Kaiser

01/2011